


How to Measure

Use a tape measure to measure around the “greater trochanter” as indicated in the diagram, then refer to the size chart. If you measure between two sizes, choose a smaller size for a firm fit and a larger size for a more generous fit.

Measure around
the Dotted line

SIZE	HIP cm	HIP ins
XXXS	79 – 86	31 – 34
XXS	86 – 91	34 – 36
XS	91 – 95	36 – 37
S	95 – 100	37 – 39
M	100 – 105	39 – 41
L	105 – 111	41 – 44
XL	111 – 118	44 – 46
XXL	118 – 123	46 – 48
XXXL	123 – 130	48 – 51



Impactactive™

HIP PROTECTORS

Slim, light and comfortable hip protectors

Disclaimer: Impactactive™ Hip Protectors are designed to help prevent a fall related hip fracture. However, while they may help lessen the severity of the impact of the fall, it does not guarantee that this will prevent a hip fracture.

Impactactive™

HIP PROTECTORS

www.impactactive.com

Order online on the website, send an email, or phone the Freephone number for assistance with your order.



info@impactactive.com

FREEPHONE

USA: 1-888-375-4890

NZ: 0800-467-228

www.impactactive.com

Impactactive™

HIP PROTECTORS

Women's Brief This ultra-comfortable feminine style has been designed with a scalloped leg trim and features an elegant embroidered waistband for that extra special touch.

Men's Boxer This classic style has been designed with a front opening and a relaxed fit leg.



Unisex Full Fit This longer leg style has been designed for those who require the extra security of a fuller fitting underwear garment. It has a fuller gusset to comfortably and securely accommodate larger pads or incontinence briefs such as pull-ups. This is also the recommended option for comfortable night-time wear.

Unisex Open Designed for those who are more frail and unable to dress or toilet themselves. The garment has velcro waistband and leg closures and opens flat to allow for easy dressing by a carer while the wearer is lying down or standing. This is also a great style for use during rehab or exercise sessions as the hip protector can be easily worn on top of clothing.

Impactactive™ and Incontinence

Both the Women's Brief and the Men's Boxer styles are designed for those who are not using incontinence products, however if a small pad is used these garments will be suitable. For those Women and Men who are incontinent - using a large pad through to a full pull-up style product, the Unisex Full Fit is the best style to choose.

Impactactive Hip Protector garments are a high quality cotton (95%) and elastane (5%) mix for a light and comfortable fit. The shields are permanently sewn into the garment and the whole garment, including the shields, can be washed and dried like normal underwear, even in commercial laundries. The shields do not need to be replaced after a fall.

HOW IT WORKS

At the moment of impact the shields momentarily stiffen.

They form a protective shield to absorb and deflect the impact force.

After impact they instantly become soft, light and flexible again.